

New equipments in the MaSC Shared Equipment Room

The following new equipments were installed.

- ① UV resin curing equipment (Room 1)
- ② Rotary polisher (Room1)
- ③ Digital camera for stereo microscope (Room3)

Please use this when sample preparation or observation is required when using shared equipment.
No other use is allowed.

Instructions for each device are located by the device.
If used, be sure to note it in the record book.

For the time being, no usage fee will be charged for the trial period.

【place】




①UV resin curing equipment (Room 1)

Put the special UV resin and the sample in the special container and turn the power on. The resin in the container will cure in about 1 minute. Convenient for processing samples.



②Rotary polisher (Room1)

Coarse to fine abrasive paper is prepared. Operation is easy with ON-OFF, direction of rotation, and two speed settings. and two speed settings.



③Digital camera for stereo microscope (Room3)

We purchased a new camera due to a breakdown of the one we had been using. It has 2 million pixels and a wide field of view and can be viewed on a PC monitor. It also has a scale gauge function. Data can be taken out by SD card.